

Foothills Prosthetics

1603 S Sterling St
Morganton, NC 28655
www.FoothillsProsthetics.com

Kevin A Wicker, CPO
kevinw@foothillsprosthetics.com
Phone: (828) 561-5164
Fax: (828) 561-5011

Ankle Foot Orthosis (AFO) Care and Use Guide

*You have just received your new orthosis which was designed and custom made especially for you.
Please refer to the following instructions for the care and use of your brace.*

Application:

Always wear a clean sock under your AFO to reduce friction and protect your skin from perspiration. It is recommended that the sock be the same length as the AFO and that you make sure the sock is wrinkle free.

Step into your AFO and tighten the straps, the step into your shoe and lace. You may need to purchase shoes that are ½ size larger to accommodate your AFO.

It is important to maintain the same shoe-heel height for which your AFO was designed. Changes in shoe-heel height while wearing your AFO may cause unneeded instability or back pain.

Wearing Schedule:

Like breaking in a new pair of shoes, it may take a brief period of time to break in your new AFO before it feels natural. To ensure a smooth process while becoming accustomed to your AFO, please follow

Skin Maintenance:

It will be necessary for you to check your skin frequently while wearing your AFO. A new AFO may cause some redness which should disappear within 15 minutes after removing the brace. If redness does not disappear, please contact Foothills Prosthetics to have this inspected.

Diabetics and those with lessened sensation are especially vulnerable to skin irritation, particularly in bony areas such as the shin or on the bottoms of the feet. Extra care should be taken, and even minor skin irritations should be treated promptly.

AFO Care:

Your AFO should be cleaned regularly with alcohol spray, or wet cloth. However, never soak the AFO as this may cause damage.



Please complete our satisfaction survey

FoothillsProsthetics.com/survey.html

If at any time your AFO becomes cracked, broken, or if you have any questions or concerns, please contact Foothills Prosthetics.

Foothills Prosthetics
Modern Technology, Old-time Hospitality®