

Foothills Prosthetics

1603 S Sterling St
Morganton, NC 28655
www.FoothillsProsthetics.com

Kevin A Wicker, CPO
kevinw@foothillsprosthetics.com
Phone: (828) 391-5164
Fax: (828) 391-5011

Care Instructions For Your Prosthesis



Silicone Liner Care:

It is highly recommended to clean your silicone liner with alcohol daily to reduce germ and residue buildup and also rotate between liners to lengthen their life span.

Daily regimen:

Turn the liner inside out and wash with warm water and a mild soap. Be sure to rinse the liner thoroughly and wipe clean with a lint free cloth. Turn liner right side out and hang or lay flat to dry.

Weekly regimen:

Turn the liner inside out and spray with rubbing alcohol. Wipe the liner clean with a lint free cloth and hang or lay flat to dry.

Prosthetic Socks Care:

It is recommended to apply clean prosthetic socks daily inside your prosthesis. Socks can be machine or hand washed with a mild soap and machine dried or allowed to hang dry.

Shrinker Care: Compression socks may be hand or machine washed with a mild soap and may be machine dried or allowed to hang dry. You will find the manufactures specific care instructions on the package.

Prosthetic Socket Care:

The inside of your prosthetic socket should be sprayed with alcohol weekly and wiped clean with a clean cloth. If at any time your socket receives any type of damage, please call your prosthetist immediately for inspection.

If your prosthesis has a locking mechanism, be sure to keep the lock clean and free of debris which could cause a lock failure. Keep your prosthetic socks in good condition; any frayed edges around the pin hole can interfere with your lock functioning properly.

Prosthetic Foot Care:

It is recommended to remove your sock and shoe weekly to inspect the prosthetic foot for any debris that may cause excessive wear.



If at any time you have questions or concerns, please contact us.

Foothills Prosthetics

Modern Technology, Old-time Hospitality®